

Bottle Compatibility Chart

Bottles						Eppendorf Rotors for Centrifuge 5920 R												Eppendorf Rotors for Centrifuge 5810/5810 R								Eppendorf Rotors for Centrifuge 5804/5804 R		
						Eppendorf Rotors for Centrifuge 5910 R																						
Bottle Manufacturer	Volume (mL)	Order No.	Size (mm) Diameter x Height	Bottom shape	Material	S-4x Universal-Large	AT*	S-4x1000 High capacity	S-4x1000 Round bucket	AT*	FA-6x250	S-4x750	AT*	S-4x Universal	AT*	S-4x500	AT*	S-4x400	AT*	S-4-104	AT*	A-4-81	AT*	A-4-62	AT*	S-4-72		
BD®/Corning®	175	352076	61 x 118	conical	PP	■	■	■	■	■	■	■	■	■	■	■	■	■**	■**	■	■	■	■	■	■	■	■	
BD®/Corning®	225	352075	61 x 137	conical	PP	■	–	■	■	■	–	■	–	■	■	–	–	–	–	–	■	–	–	–	■	–	■	
Nalgene®	175	3143-0175	61.5 x 143.5	conical	PP	■	–	■	■	■	–	■	–	■	–	–	–	■***	–	■	–	–	–	–	–	–	–	
Nalgene®	175	3144-0175	61.5 x 144.5	conical	PC + Screw closure	■	–	■	■	■	–	■	–	■	–	–	–	■***	–	■	–	–	–	–	–	–	–	
Nalgene®	175	3145-0175	61.5 x 135.6	conical	PS + sterile	■	–	■	■	■	–	■	–	■	■	–	–	■***	■***	■	–	–	–	–	–	–	–	
Nalgene®	250	3140-0250	61.5 x 135.1	flat	PC	■	–	■	■	■	■***	■	–	■	■	–	–	■**	■**	■	–	–	–	–	–	–	–	
Nalgene®	250	3127-0250	60 x 128.8	flat	FEP	■	–	■	■	■	■	■	–	■	■	■	–	■**	■**	■	–	■	–	■	–	■	■	
Nalgene®	250	3122-0250/ 3120-0250	61.8 x 127.6	flat	PP	■	–	■	■	■	■	■	–	■	■	■	–	■**	■**	■	–	■	–	■	–	■	■	
Nalgene®	500	3141-0500	69.5 x 160	flat	PP + sealing cap	■	–	–	■	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	
Nalgene®	1000	3122-1010	98.1 x 180.2	flat	PC	–	–	–	■	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	
Nalgene®	1000	3120-1010	98 x 179	flat	PP	–	–	–	■	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	
Nunc®	200	376813	60 x 130	conical	PP	■	–	■	■	■	–	■	–	■	■	■	–	■***	■***	■	–	■	–	■	–	■	–	–
Nunc®	250	376814	60.2 x 145	conical	PP	■	–	■	■	■	–	■	–	–	–	–	–	–	–	■	–	–	–	–	–	–	–	
Corning®	250	430776	60 x 164	conical	PP	■	–	–	■	■	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	
Corning®	250	431841	61 x 126	flat	PP	■	■	■	■	■	■	■	–	■	■	■	–	■**	■**	■	–	■	–	■	–	■	■	■
Corning®	250	431842	61 x 128	flat	PC	■	–	■	■	■	■	■	–	■	■	■	–	■**	■**	■	–	■	–	■	–	■	–	■
Eppendorf®	400	5810 729.009	80 x 133	flat	PP	–	–	–	–	–	–	–	–	–	–	■	–	–	–	–	–	–	■	–	–	–	–	
Eppendorf®	500	5820 707.003	83 x 134	flat	PP	–	–	–	–	–	–	–	–	–	–	■	■	–	–	–	–	–	■	■	–	–	–	
Eppendorf®	750	5820 708.000	102 x 133	flat	PP	■	■	–	■	■	–	■	■	■	–	–	–	–	–	–	■	–	–	–	–	–	–	
Eppendorf®	400	5910 707 002	81 x 117	flat	PP	–	–	–	–	–	–	–	–	–	–	–	–	■****	■****	–	–	–	–	–	–	–	–	–
Eppendorf®	1000	5910 770 006	118 x 133,1	flat	PP	–	–	–	–	–	–	–	–	■****	■****	–	–	–	–	–	–	–	–	–	–	–	–	–
Corning®	500	431123	95 x 149	conical	PP	■	■	–	■	■	–	■	–	–	–	–	–	–	–	–	■	–	–	–	–	–	–	
Corning®	500	431844	–	flat	PP	■	–	–	■	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	
Corning®	500	431847	–	flat	PP	–	–	–	■	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	
Beckman Coulter®	500	355605	69 x 160	flat	PC	■	–	–	■	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	
Beckman Coulter®	500	355607	69 x 160	flat	PP	■	–	–	■	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	
TPP®	600	87600	100 x 183	conical	PP	–	–	–	■	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	–	

Technical specifications subject to change. Errors and omissions excepted.

■ Yes * Aerosol tight cap possible ** Use with adapter 250 mL flat bottom shape only. *** Use with adapter 175-225 mL conical bottom shape only. **** Available soon

www.eppendorf.com/centrifugation

Your local distributor: www.eppendorf.com/contact
 Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany
 E-mail: eppendorf@eppendorf.com · www.eppendorf.com

Beckman Coulter® is a registered trademark of Beckman Coulter, Inc., USA. Corning® is a registered trademark of Corning, Inc., USA. Nalgene® is a registered trademark of Nalge Nunc International Corporation. Nunc® is a registered trademark of Nunc A/S Corporation, Denmark. TPP® is a registered trademark of TPP Techno Plastic Products AG, Switzerland. Eppendorf® and the Eppendorf Brand Design are registered trademarks of Eppendorf AG, Germany. All rights reserved, including graphics and images. Copyright © 2019 by Eppendorf AG.